

# UNCLASSIFIED

AD NUMBER
ADB229914
NEW LIMITATION CHANGE
TO Approved for public release, distribution unlimited
FROM Distribution authorized to U.S. Gov't. agencies only; Specific Authority; 3 Oct 97. Other requests shall be referred to U.S. Army Medical Research and Materiel Command, Fort Detrick, MD 21702-5012.
AUTHORITY
USAMRMC ltr, 4 Dec 2002

THIS PAGE IS UNCLASSIFIED

AD \_\_\_\_\_

CONTRACT NUMBER DAMD17-97-C-7027

TITLE: Rapid Membrane Based Immunoassays for Detection of Toxins and Infectious Disease Pathogens

PRINCIPAL INVESTIGATOR: John Colanduoni, Ph.D.

CONTRACTING ORGANIZATION: Arista Biologicals, Incorporated  
Bethlehem, Pennsylvania 90501

REPORT DATE: July 1997

TYPE OF REPORT: Final, Phase I

19971002 028

PREPARED FOR: Commander  
U.S. Army Medical Research and Materiel Command  
Fort Detrick, Maryland 21702-5012

3 Oct 97

DISTRIBUTION STATEMENT: Distribution authorized to U.S. Government agencies only (specific authority). Other requests for this document shall be referred to U.S. Army Medical Research and Materiel Command, 504 Scott Street, Fort Detrick, Maryland 21702-5012.

The views, opinions and/or findings contained in this report are those of the author(s) and should not be construed as an official Department of the Army position, policy or decision unless so designated by other documentation.

DTIC QUALITY INSPECTED 4

# REPORT DOCUMENTATION PAGE

Form Approved  
OMB No. 0704-0188

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington, VA 22202-4302, and to the Office of Management and Budget, Paperwork Reduction Project (0704-0188), Washington, DC 20503.

1. AGENCY USE ONLY (Leave blank)		2. REPORT DATE July 1997		3. REPORT TYPE AND DATES COVERED Final, Phase I (25 Nov 96 - 24 May 97)	
4. TITLE AND SUBTITLE Rapid Membrane Based Immunoassays for Detection of Toxins and Infectious Disease Pathogens				5. FUNDING NUMBERS DAMD17-97-C-7027	
6. AUTHOR(S) Colanduoni, John, Ph.D.					
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Arista Biologicals, Incorporated Bethlehem, Pennsylvania 90501				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING / MONITORING AGENCY NAME(S) AND ADDRESS(ES) U.S. Army Medical Research and Materiel Command Fort Detrick, Maryland 21702-5012				10. SPONSORING / MONITORING AGENCY REPORT NUMBER	
11. SUPPLEMENTARY NOTES					
12a. DISTRIBUTION / AVAILABILITY STATEMENT Distribution authorized to U.S. Government agencies only (specific authority). Other requests for this document shall be referred to U.S. Army Medical Research and Materiel Command, Fort Detrick, Maryland 21702-5012.				12b. DISTRIBUTION CODE 3 Oct 97	
13. ABSTRACT (Maximum 200 words)  Arista Biologicals, Inc. proposes to develop a series of assays to detect infectious disease pathogens and toxins on a series of rapid nitrocellulose membrane based assays utilizing colloidal gold and fluorescent dyes in both a unitized and multi-analyte format capable of semi-quantitation with a simple, inexpensive reader similar to a home diabetic glucose meter. We have formatted and initially tested 6 of the pathogens in the unitized format and have developed working tests for both dengue and hanta viruses. We plan to continue our investigation of the use of fluorescent dyes side by side with visible dyes, as the initial work has shown as at least comparable results and indicates that with a UV reader we may well find the enhanced sensitivity that we originally predicted. We have completed the exploratory work and initial testing of a reader system for enhanced readability and quantification and plan on continuing those efforts in the Phase II segment of the project. All of these efforts will enhance the detection capabilities of rapid assays for both military forward deployment and civilian point of care applications					
14. SUBJECT TERMS Immunoassay, Quantifying Reader System, Environmental Testing, Toxin Detection, Disease Detection				15. NUMBER OF PAGES 65	
				16. PRICE CODE	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT Limited		

## FOREWORD

Opinions, interpretations, conclusions and recommendations are those of the author and are not necessarily endorsed by the U.S. Army.

( ) Where copyrighted material is quoted, permission has been obtained to use such material.


( ) Where material from documents designated for limited distribution is quoted, permission has been obtained to use the material.


( ) Citations of commercial organizations and trade names in this report do not constitute an official Department of the Army endorsement or approval of the products or services of these organizations.

( ) In conducting research using animals, the investigator(s) adhered to the 'Guide for the Care and Use of Laboratory Animals,' prepared by the Committee on Care and Use of Laboratory Animals of the Institute of Laboratory Animal Resources, National Research Council (NIH Publication No. 86-23, Revised 1985).

( ) For the protection of human subjects, the investigator(s) have adhered to policies of applicable Federal Law 32 CFR 219 and 45 CFR 46.

( ) In conducting research utilizing recombinant DNA technology, the investigator(s) adhered to current guidelines promulgated by the National Institutes of Health.

  
Principal Investigator's Signature

  
Date

## **TABLE OF CONTENTS**

<b>FRONT COVER</b>	<b>P.1</b>
<b>REPORT DOCUMENTATION PAGE (SF 298)</b>	<b>P.2</b>
<b>FOREWORD</b>	<b>P.3</b>
<b>INTRODUCTION</b>	<b>P.5</b>
<b>THE PROJECT</b>	<b>P.7</b>
<b>GENERAL METHODOLOGY</b>	<b>P.7</b>
<b>GENERAL INFORMATION</b>	<b>P.10</b>
<b>PROJECT PROGRESS REPORT</b>	<b>P.11</b>
• <b>FLUORESCENT PROBES</b>	<b>P.11</b>
• <b>HAND-HELD READER SYSTEM</b>	<b>P.13</b>
• <b>MULTI-ANALYTE FORMAT</b>	<b>P.14</b>
• <b>HANTA VIRUS &amp; DENGUE VIRUS ASSAYS</b>	<b>P.21</b>
<b>CONCLUSION</b>	<b>P.24</b>

### **ATTACHMENTS:**

- **ATTACHMENT A - SAMPLES OF READER GRAPHING**
- **ATTACHMENT B - REQUEST TO AMEND CONTRACT**
- **ATTACHMENT C - SAMPLE PACKAGE INSERT**
- **ATTACHMENT D - INFORMATION ON (+) DENGUE SERA**
- **ATTACHMENT E - METHODOLOGY PROCEDURES**

## **INTRODUCTION**

Arista Biologicals, Inc. (ABI), as the prime contractor, and Majesco Biologicals, Inc., as an affiliated company, propose to develop a series of rapid membrane based tests to screen medical and environmental samples for the presence of toxins, bacterial antigens, and/or antibodies and to develop a simple reader system capable of increasing sensitivity and specificity, enabling quantification, and removing objectivity in visual reading of such membrane based tests. The aim is to address the problems experienced with diagnostic immunoassay testing in the field and/or by laymen and less technically personnel.

There has been a stereotyping of these tests as being somewhat inferior to more classic methods of diagnostic testing. Many early generation tests were inferior in quality and prematurely released before proper controls were available. This led many in the medical field to believe these tests were not capable of the same accuracy as ELISA, PCR or biosensors. It is our belief that by using the basic immunochromatographic membrane based assays it is possible to develop tests equal, or superior, in accuracy and dependability as these other types of assays. Most of the other confirmatory methods require more expensive equipment, sophisticated electronic technology and long multi-step procedures, drawbacks not associated with our assay format.

These benefits are important features for a test that is intended to be performed by users in the field with limited technical training. The system's design will make it much less expensive, much faster and easier to perform by the end user than other technologies including ELISA, PCR and bio-sensor systems. In that ELISA is the test most frequently compared to our rapid membrane format, it should be pointed out that the differences are significant:

### **ELISA**

Multi-step  
Minimum of 30-60 mins.  
Primarily instrument based  
labile stability  
cost effective if multiple specimens  
much higher tech know-how  
single analyte

### **Membrane**

One step  
3-5 minutes  
non-instrumented  
totally stable for 2yrs  
same cost single or multiple  
anyone can run  
potential for multiple analytes

Currently, ABI utilizes colloidal gold particles as the conjugate in the above described chromatographic system. In the past, we have successfully produced a variety of assays for infectious disease testing. As part of our Phase I project, we investigated the use of fluorescent probes to replace or enhance the colloidal gold system, either attempting to attach the fluorescent probes directly to the detector antibody or proteins and/or attaching the fluorescent probes to the colloidal gold detector protein conjugates. The use of any such fluorescent probe modification should enable us to increase the sensitivity of these membrane based assays as much as one thousand fold. This adaptation can increase the

use of this technology not only as a screening test but as a confirmatory assay and will allow that the system will still be visually readable in situations when the optional reader with its enhanced sensitivity and objectivity features are not required.

As a second goal of the Phase I SBIR, ABI has modified this immunochromatographic assay into a multi-analyte test format that will make it more applicable in both the military and civilian sectors such as in detection of environmental contaminants, and screening to differentially diagnose infectious diseases.

The potential benefits associated with these assays and reader systems have garnered the interest of both environmental and infectious disease sections of the DOD as well as that of several well placed commercial partners. These assays have particular application for military forward deployment and point of care use in Third World countries where the assay characteristics of simple, one step, non-instrumented, non-technical expertise, 3-5 minute result assays are essential. This presents a creditable scenario for the Dual Use role of these assays. The exceptional stability of the components/reagents of these products under any conditions is an additional factor in their desirability in all markets.

Arista staff have developed assays, in house, for individual clients whose interests cover most areas of diagnostics, such as veterinarian, environmental, drug testing, bacterial and viral pathogens, hormones, specific antibody, etc.. ABI staff also train other companies in the production of these types of tests for either standard products or products developed specifically for the client. Our staff has up to ten years experience in this area of the business. The assays developed can be used for serum, plasma, whole blood, urine, feces, airborne particles, dirt, etc. Several of the assays have been tested in trials, including very successful field trials of plague, anthrax, and leishmania.

## **THE PROJECT:**

This project has three basic areas that are being investigated simultaneously:

1. The first task is to develop a rapid, hand-held immunoassay system capable of being used in the field with state of the art sensitivity utilizing fluorescent probes combines with colored particles in an effort to enhance detection capability and allow visualization by an ultraviolet light source.
2. The second task is to develop an immunoassay system capable of being read with a small hand-held reader that will perform quantification, providing a permanent written record. This system will also allow for visual determination of sensitivity utilizing a densitometer.
3. The third task in this project will be the ability to perform several assays on a single device in order to reduce the number of individual tests being run by the operator. These multi-analyte tests would differentiate between multiple agents that are suspected of being present.
4. As a component of this project, work has been done to optimize the basic format for the detection of Hanta virus and Dengue virus with the ability to develop a test for relevant antibody detection.

## **GENERAL METHODOLOGY:**

In this section we provide some detail on specific materials/ reagents used within our tests as well as some narrative on methodology for preparing actual test strips in our format. A more complete description of production procedures and methodology is contained in a separate, 40 page document submitted as an attachment to this report. Note that unless there is some variation in procedures, methods, or reagents, this information will be presented only in this section and not for each of the individual projects.

- **buffer optimization** - Assay buffer is incorporated into the sample pad. We start with standard pad for serum and buffer and work out from there. This pad is soaked with 200 mM Sodium Borate pH 9.3. 0.1% Triton X-100, 0.5% PVA. We will normally adjust this formulation based on whether we see false positive results. For example, we will raise the pH up to as high as 11 and increase the Triton X-100 to up to 1% or add up to 1% SDS. If false positives are still a problem we will add an assay buffer such as 1% BSA, 0.1% Sodium Azide, pH 7.4 50 mM Sodium Phosphate and dilute the sample with that diluent directly on the assay pad, thereby diluting out the serum as much as




75%. We can also add any specific reagent to remove problems such as E. coli extract, goat serum, other surfactants, etc. This is normally only necessary for antibody assays where the antigen we spray on the membrane contains contaminants which the sera sometimes reacts with. The above described diluent may also be used as a general diluent to increase flow rate of viscous or lipemic samples which when run neat sometimes take over 10 minutes to run.

- **membranes (pore size) and colloidal gold (particle sizes)** - Our previous experiences with other strip colloidal gold systems have taught us that a pore size of approximately 15 microns is the best for both flow rate and reactivity. If larger pore sizes are used the flow rate is usually too fast for reactions to take place at low sensitivities and if smaller pore sizes are used a test is difficult to finish in 5 minutes. We generally use the SR (15 $\mu$ ) series of nitrocellulose membranes from Millipore. We did, however, look at 20 micron pore size nitrocellulose for the fluorescent latex assays used in this study because the latex is about 10 times larger than the colloidal gold. As to colloidal gold particle size, previous experience in gold assays has shown us that a 40 nm particle is the best for the strip assays. Smaller particles give less signal because of how the gold scatters light and larger particles tend to migrate slower, generate a purple blue color, and are difficult to work with (stability). The latex had to be approximately 250 nm in order to be able to centrifuge and process the conjugates. This is due to its density being approximately 9 times less than that of the gold.
- **production methodology** - The first step is to take 30 cm plastic cards (3 inches by 30 cm) and peel off the cover of the double stick adhesive that is present 1 inch lengthwise in the center of the card. The nitrocellulose membrane is then laminated onto this section prior to dispensing the capture reagent. The capture reagent is then dispensed on the card using the Biodot Quantijet Dispenser at a volume of approximately 1 microliter per cm at a predetermined concentration (typically at 1 to 4 milligrams per ml). The card is dried at 37 degrees for one hour, then normally blocked with a BSA sucrose (20% sucrose) mixture, and then re-dried as before. While the cards are drying the conjugate pad material, polyester, is also pretreated with PVA and surfactant to make it hydrophilic and allow conjugate to dry and later re-hydrate. The absorbent pad of cellulose acetate is also soaked with the assay buffer (normally made up of sodium borate, BSA and or PVA, and surfactant(s)). These materials are also air dried and can be stored for long periods of time. The colloidal gold or latex conjugate is then dispensed on the conjugate pad material at the predetermined concentration and volume and is then also air dried at 37 degrees Celsius. After all of the materials are dried, the absorbent pads, conjugate pads, and sprayed and blocked cards are assembled by removing the rest of the covers of the double stick tape on each card (overlapping the conjugate pad on the bottom of the nitrocellulose and the absorbent pad which itself overlaps with the conjugate pad). Untreated

cellulose acetate is then placed and overlapped with the upper portion of the nitrocellulose. The cards are then slit into strips of the desired width, (this width is 4 cm for the device used in this project) and stored with desiccant while awaiting final assembly. The strips are later placed into the plastic housings and closed and pouched into foil pouches along with desiccant and droppers. Pouches are labeled with lot # and name of test. Sensitivity and specificity is tested for each lot to assure appropriate antigen or antibody levels of detection.

(Figure 1)

Standard test format



absorbent pad - cellulose acetate  
sample pad - cellulose acetate

conjugate pad - polyester  
conjugate - protein A colloidal gold (PAG)

- **format methodology** - the sample is added to the sample well (**Figure 1**):
  - The sample reconstitutes the assay buffer. The filter materials also acts to filter out any large particles that might contaminate the sample. The next material that the sample comes in contact with is the conjugate pad (containing the identifying antibody). The conjugate pad is also buffer treated so as to promote reconstitution of the conjugate and maintain immunological activity of the dyed or colored conjugate.
  - The test sample then can react with the conjugate if the test sample contains antigen specific for the antibody conjugate.
  - This mixture then migrates onto the nitrocellulose membrane. The nitrocellulose membrane section has two horizontal lines of reagent material sprayed onto it in a single analyte format and multiple horizontal lines in the multi-analyte format.

- When a capture protein or related material is sprayed onto the nitrocellulose, it becomes non-covalently attached to the membrane because of the membrane's high affinity for proteins and other hydrophobic characteristic materials. The lower band (or bands in multi-analyte) has (have) the capture antibody sprayed onto it in a sandwich assay or the specific antigen sprayed onto it in an antibody capture assay. The conjugate has the detector antibody attached to it in a sandwich assay and the conjugate has the "anti-antibody" reagent attached to it in the antibody capture assay. If the conjugate has the correct antigen attached to it so that the capture protein recognizes it, the colored or dyed particle will then attach to the lower line of the strip creating a visual or visualizable horizontal line.
- Any conjugate that migrates beyond the first horizontal test band will then attach to the second band since the second band is always an anti-conjugate material. For example, in the case of a monoclonal antibody conjugate an anti-mouse IgG antibody is used as a positive control capture reagent. Therefore, two horizontal bands are seen for a positive reaction and only one, the higher band, is observed in a negative test result.
- If for some reason the test is invalid, there will be no band observed and the sample will need to be re-tested or investigated for the cause of the invalid result.

#### **GENERAL INFORMATION:**

- **separation of IgM and IgG** - We have performed preliminary experiments to separate out reactions to IgM and IgG simultaneously on a single strip for the antibody assays. We decided to use the test for dengue as the prototype. We did this by placing two test bands on the strip. One was goat anti-human IgM and the other anti-human IgG. Dengue antigen was conjugated to the colloidal gold particle. Previous experience with Epstein-Barr assays had shown us that this is difficult to accomplish with good success, however in that assay we were working with lysate type antigen and we thought that the recombinant antigen would work better in this manner. We have seen some separation of signal, but it definitely needs more work, which we will continue in Phase II.
- **avidin/biotin use** - A streptavidin/biotin common reagent system was evaluated but was not found to enhance sensitivity and was not even as sensitive as direct conjugation of primary antibodies. We had similar experiences earlier on in other projects. This is probably because our assay is very similar to a homogeneous assay and is not sequential in nature where streptavidin biotin works best. We have also had problems with this system

with individuals who were on large doses of vitamins where it interfered with biotin binding because of similarities in structure with biotin of some vitamins.

- **use of F-1 and plague antigens** - We are preparing to immunize with plague antigens. We would need to commit additional goats but do not see spending money till we see if we are going to continue, in that there is limited use for our own commercial applications. We also need more antigen for making affinity columns. Antisera will be useless without antigen for affinity purification. Antibody material we have from the Army now is good enough, as far as activity levels, for a feasibility study.

## **PROJECT PROGRESS REPORT:**

### **1. Fluorescent Probes**

**Strategy 1** - The fluorescent probes were evaluated using several methods in our attempt to increase sensitivity of the strip assays. Because of the nature of the experiments these trials used hundreds of milligrams of antibodies. The target reagents (ricin, SEB, etc.) are available in only limited quantities, and therefore hCG reagent (in abundant, inexpensive supply) was used as the model. Initial studies used direct conjugation of the fluorescent tags to the monoclonal antibody replacing the colloidal gold conjugates normally used in the strips. Rhodamine, Texas red, and fluorescein were used as the fluorescent conjugate. The conjugation **methodology** was as follows:

The initial experiments conjugating the fluorescent probes directly to the antibodies were all performed using isothiocyanate derivatives of the organic compounds. Isothiocyanate will couple directly to proteins without requiring other chemicals. We conjugated in pH 7.4 PBS at 10 mg of probe to 1 mg of protein. They were allowed to stir at 4 degrees centigrade in the dark overnight and dialyzed three times at 1000 fold excess dialysis buffer. The dialysis buffer was pH 7.4 PBS.

The highest color intensity visually observed when ultraviolet light was placed on the strips was observed when using fluorescein. Fluorescein was attached using FITC as the fluorescent tag and optimal conjugation observed at approximately 10 mg of FITC to 1 mg of antibody at pH 7.4 in PBS. However, when placed in the standard hCG strip it appeared to be only as sensitive as the normal colloidal gold assay to approximately 1 ng/ml. We also needed to use approximately ten times more antibody conjugate per test using about 10 ug of antibody per test compared to a colloidal gold antibody conjugate.

**Strategy 2** - The next strategy taken was to add the fluorescent probe to the colloidal gold antibody conjugate directly hoping that we would get a higher FITC to particle ratio than conjugating to individual antibodies since the antibody gold conjugate contains many

antibody molecules per particle. The gold is approximately 40 nm in size and coated with many smaller antibody molecules. When tested, this showed almost no fluorescence observed at even 100 ng/ml levels of hCG. Several levels of hapten conjugation were attempted from 1 to 100 mg of label per mg of conjugate without any positive results. This may be due to interactions with the more hydrophobic gold or contaminants in the gold preparations.

**Strategy 3** - The final, and most successful method, was found to be with using latex fluorescent beads and coating the hCG antibody to these particles in the following manner. Three fluorescent types of beads were purchased with approximately 300 nm particle size. The 300 nm was chosen because at smaller size particle ranges it is very hard to separate unconjugated antibody from the antibody bead molecules since the density of the latex particles is not much higher than water. This compares to the gold particles which have an approximate tenfold higher density than water. Larger latex particles tend to clog the nitrocellulose membranes being used in this project since they are close to the membrane pore sizes used (typically 3 to 15 microns) and also tend to settle in solution when being applied.

Fluorescein, Texas Red, and Rhodamine containing latex particles were purchased. Both standard polystyrene and carboxylated particles were purchased for antibody conjugation optimization. Conjugation **methodology** for the latex particles was as follows:

Using approximately 300 nm polystyrene latex particles from Bangs Laboratories and using FITC, Texas Red, rhodamine dyed latex, the latex particles were coated in the following manner. The latex particles were diluted to 0.1% solids/volume in pH 6.8 100 mM Sodium Phosphate. Antibody was added at 1 mg/10 ml of diluted latex particles. The mixture was rocked overnight at room temperature. The latex was washed three times by centrifugation and re-suspension into pH 7.4 PBS in a tabletop centrifuge at 3000 rpm for 15 minutes and then re-suspended by vortex. The final concentration was 1 ml per 10 ml of starting material making the final concentration 1% solid/volume. The conjugates were used directly in the assays in both a liquid and dried form on the strip.

Latex particles carboxylated were conjugated in the following manner. The latex was diluted into PBS to 0.1% solid/volume and solid carbodiimide added to the mixture and rocked for 30 minutes for activation. The mixture was then centrifuged for 15 minutes at 3000 rpm and re-suspended twice into pH 7.4 PBS. Antibody was added as above and washed as above. The latex thus produced seemed identical to that in the first procedure so that the final decision was to use the polystyrene latex and the passive absorption method. We had hoped the 2<sup>nd</sup> method would allow for an increase of binding or more activity due to covalent attachment instead of passive absorption.

In this latex phase of the project, it was found that again the fluorescein type materials showed the highest fluorescence. We expected the carboxylated particles to work better but, suprisingly, this was not the case. Normally, covalent attachment forms a non-

reversible binding making a more stable conjugate and a more active conjugate than passively absorbed conjugates using the nude polystyrene particles. This is thought to occur since the covalent attachment allows for more free movement of the antibody on the surface of the particle. The conjugates were then placed on strips coated previously with hCG, varying the hCG in the sample and the volume of conjugate. The passively absorbed material was found to be about 5 times more active. The sensitivity again was found to be approximately 1 ng/ml although now we needed approximately 1 microgram of antibody per test.

Although initially we were disappointed to see that we did not detect the expected increase in sensitivity using these fluorescent probes as compared to the colloidal gold, discussing it with several peers in the industry allowed us to arrive at the **following conclusion**. Normally, when using fluorescent probes in immunoassays, the fluorescence is not visually read but rather is read by instrumentation. The signal may be in the non-visual range towards the ultraviolet range of the spectrum. We most likely need to use instrumentation as discussed in the original proposal. The type of reader needed is more complicated than the visual readers studied during this phase I part of the project and we are currently exploring the possibility of prototyping such an instrument with the appropriate engineering experts.

**NB:** Diagrams of format are identical to gold assay, with substitution of fluorescent conjugate for gold. Methodology for making the test is also identical to gold assay. Therefore, we will not repeat the information or the diagrams.

## **2. Hand-Held Reader System**

As noted in our prior reports, we have been working with a reader obtained from Source Instruments. This reader was used in the following manner:

The strip was placed in the tray at a preset position (for reproducibility). This meter scans a 1 cm portion of the strip. The data recording procedure provided was loaded onto a PC. The data recorded are position versus intensity measurements. These, as well as any other parameters for data intake, are not adjustable with this specific instrument. For example, the zero is predetermined internally, and we cannot offset background colors or increase sensitivity. This is the major drawback with this prototype since we cannot alter data input at all. The strip is then scanned automatically when the start button is hit. The data then is stored in the PC for later manipulation. We are thus provided with the computer graph intensity versus length and must then integrate the area under the curve for the test line. We need to use our own program for this as the manufacturer does not have integration software available at this time. After collecting data at several antigen levels and doing the data integration, we then plot the integrated signal versus antigen concentration. The plot shows sigmoidal type results at the low end of the curve. This is probably due to the signals' relative intensity being measured. **(See Attachment A for sample graphing)** We probably should take the log function and plot that so as to be measuring absorbency. However, since we cannot zero the instrument, we decided not to

attempt that procedure at this time.

We have identified and obtained an additional prototype color reader for evaluation. This reader comes from Biotronics Associates. It differs from the Source reader in several important aspects. It is adjustable for negative readings and it can be read both by absorption or transmission. The fact that you can re-zero the reader to different negative backgrounds makes it more reliable when the background of the strips are different such as when using serum samples with different background color or using strips with different gold conjugate levels (which may cause a higher or lower reddish color background across the entire strip). The wavelength, using this new reader, is also adjustable so that you are not looking at a wide range of colors. The earlier Source reader was not color adjustable and had an arbitrary zero level internally set which was not adjustable. Both the Biotronics and the Source readers are sensitive to the nanogram range using the hCG assay as a standard. However, the color intensity was linear with the Biotronics Associates reader in the entire range. The Source reader was observed to be nonlinear in the lower ranges. The Biotronics reader also has a direct readout and does not require a PC. The Biotronics reader does not scan but has a 5 mm orifice within which the line on the strip can be placed. The zero can be calibrated to a blank spot on a strip to correct for background color differences. We can get linear results measuring absorbency but find that great care must be taken in placing the strip correctly over the light orifice. A scanning method would seem preferable with the strips but using a more narrow orifice might alleviate some of the problem with the Biotronics. The line is approximately 1 mm in thickness, at this time, and we are using a 5 mm wide strip for these experiments.

We feel that we can combine the positive aspects of both readers and incorporate them into the next generation of instruments. Both instruments are approximately 10 inches square in shape 4 inches in height. The final prototype needs to approach hand-held size. Our experience with manufacturers that we have worked with makes us confident in the ability to progress further on an appropriate prototype.

### **3. Multi-analyte format**

Our initial proposal included a plan to format the multi-analyte in a radial design. This, we felt, had potential for easy workability and readability. (For edification of designs, please see **Figures 2-7**). We optimized the radial assay format using both hCG and drugs of abuse assays. We used hCG and drugs of abuse because of their relevant ease of being obtained. The high volumes of reagents necessary for this procedure relegate the use of target reagents (SEB, ricin, plague, VEE, etc.) to only the final stages of development. We ran into significant problems with the radial assay which, at this point in time, have alerted us to the need to abandon it, if only temporarily, in favor of a strip format. It is certainly a possibility to go back to re-examine the radial design as part of a Phase II program. The drawbacks that we ran into with the radial format were:


- we required approximately 25 fold more conjugate per assay using this radial format

since most of the conjugate does not migrate specifically to the target assay spot.

- The sensitivity of the assays was diminished when we used a common anti-mouse IgG conjugate with free individual monoclonal antibodies for each specificity instead of a different conjugate for every assay. Therefore we needed to use an excess of conjugate for every assay using conjugate of each specificity.
- We found that we used approximately 30 times the amount of the nitrocellulose membrane which is a relatively expensive item. The circular/radial nature of this format requires us to die cut out circles from sheets (rather than the narrow strips normally used where rolled stock can be used directly). This wastes even more membrane and consumes production time.
- The cost per test for the membrane alone is approximately 90 cents per test, versus 40 cents in the strip format.
- We also found that the reagent dispensing and processing equipment we currently have is not sophisticated enough to use directly in the radial design and that much more sophisticated production equipment would be necessary. The approximate cost for such equipment is estimated to be \$750,000 and would be required to be custom made.



**Figure 2**      **ANTIGEN DETECTION**  
**RADIAL ASSAY**



**Figure 3**      **ANTIBODY DETECTION**  
**RADIAL ASSAY**




Figure 4

# ANTIBODY ANTIGEN ANTIBODY GOLD COMPLEXATION

## EXAMPLE 1

antigen 1 and 2 positive  
antigen 3 negative





Figure 5

# SAMPLE RECONSTITUTED CONJUGATE

## EXAMPLE 1

antigen 1 and 2 positive  
antigen 3 negative





Figure 6

# ANTIGEN ANTIBODY PROTEIN A COMPLEXATION EXAMPLE 2

antibody 1 and 2 positive  
antibody 3 negative







Figure 7

# SAMPLE RECONSTITUTED CONJUGATE EXAMPLE 2

antibody 1 and 2 positive  
antibody 3 negative




### Multi-analyte format (Con't)

For all of these reasons we decided to explore the multi-analyte strip format for the remainder of the project. This format may be more useful for antibody assays than antigen assays, in the future, since only a single conjugate would be used for all the specificities on the device. As a first step in this multi-analyte format design, we optimized, in individual strip format, tests for ricin, SEB, plague, and VEE. **(Figure 8)**

**(Figure 8)**

#### VEE Assay \*\*



**\*\*Exact same format for ricin, plague, and SEB**

The sensitivities of these current generation of the individual assays were acceptable, although we would like to work to increase sensitivity.

The next step in the multi-analyte strip development was to spray all four assays on an individual strip and optimize the mixture. **(Figure 9)** We have observed some nonspecific binding of one colloidal gold conjugate with another capture antibody and need to optimize conditions where these antibodies and conjugates do not have any interactions. We see this problem frequently with assays made from mixtures of antibody from different source animals where there is heterophile type interactions in which one species has antibodies or specific interactions with another species. Normally each assay is optimized individually and the buffer in the absorbent pad and blocking of membranes for each assay are optimized in separate procedures to prevent these type of interactions.


Other, more tedious methods, for reducing this type of nonspecific interaction are:

- to individually absorb out these heterophile interactions using agarose antibody columns to remove these antibodies from the already affinity purified polyclonal antibodies.
- to cleave either the monoclonal or the polyclonal antibodies and form Fab type fragments since removing the Fc region of the antibodies tends to remove these interactions also.

We are in the process of first optimizing the buffers before we use either of these two modifications of the antibody treatments since the amount of reagent available is relatively small. Because of scaling-up problems associated with these procedures, both these protein purification and modification methods tend to be reagent consuming and are only feasible with quantities over 100 mg of purified reagent.

An example of a multi-analyte format for antibody detection (using hanta and dengue antigens) is provided in **Figure 10**.

(Figure 9)



**KEY:** A = mixture of monoclonal anti - VEE, ricin, SEB, plague colloidal gold conjugates

B = polyclonal anti-plague band


C = polyclonal anti-SEB band

D = polyclonal anti-ricin band

E = polyclonal anti-VEE band

F = anti IgG antibodies band (control)

(Figure 10)



#### **4. Hanta virus and dengue virus assays**


##### **4a. Hanta virus assay (Figure 11)**

In our efforts to develop a prototype for this test, we have evaluated two sets of hanta antigens furnished by the Army and have found, in our system, that they are immunologically inactive, testing against currently available clinical samples of positive sera obtained from USAMRIID. We have investigated the possibility that there are chemicals in the system which would make them incompatible with our nitrocellulose binding and have concluded that this is not the problem. These antigens must lose their immunological activity upon binding to the membrane, be contaminated with an inactive protein or other form of the antigen, or have very low immunological activity relative to antigens used successfully in this type of assay format. We are currently waiting for additional material from the Army.

We have found, however, that a recombinant antigen, prepared in *E. coli*, from a private South Korean company (Che Dang, Seoul, S. Korea) seems to work very well. When using this recombinant antigen we can produce a test which gives us excellent specificity and sensitivity. By coating approximately 500 ng of this antigen per test strip and using a protein A colloidal gold we have been able to consistently detect antibody positive test samples. Clinical samples available from domestic sources are relatively rare although we have observed positive results with the two samples made available to us by USAMRIID. The only information we have on those 2 samples is that they should be strongly positive but are quite old. Ten problematic type negative samples used in-house have shown negative results (i.e., no false positives). Since they have access to many more samples at the Korean NIH than we do here in the States, we have sent 90 tests to Korea for their evaluation. They have found that out of the 70 positive samples tested thus far they observed 68 as clear positives and had not observed any false positives.

(Figure 11)

### Hanta Virus Assay



In the interim, we have made several improvements to increase sensitivity and have provided Korea with approximately 120 more assays for their continued evaluation. We have also supplied the Army with samples of both the raw Korean antigen and tests for evaluation at USAMRIID. The antigen supplied to us by the S. Koreans is estimated to be only about 50% pure. They are preparing a more pure antigen for us, to see if there is any further improvement that can be achieved in the assay.


**Specific methodology for hanta assay** - The hanta test was formatted in the following manner. The antigen was coated at 1 mg/ml at 1 microliter per cm on SR Millipore nitrocellulose membrane in PBS. The card was blocked with BSA Sucrose PBS and air dried. The colloidal gold used was a Protein A conjugate with a 40 nm particle size at an OD of 50 at 2 ul/cm. This was the maximum concentration of antigen available at this time. The assay buffer was applied to the absorbent sample pad. The buffer found optimally to work was pH 9.3 200 mM Sodium Borate, 0.1% BSA, 1% Triton. The assay was found to require the addition of E coli extract (provided from the Korean antigen source) at 10% in passive gold diluent. This mixture was found to be necessary in at least a 25% ratio to the total volume of the sample in order to remove all false positives from problem negative samples. Evidently enough E. coli debris is left in this antigen preparation to give false positives with individuals with a titer to E coli.

#### **4b. Dengue Assay (Figure 12)**

As part of our attempt to develop a working assay for dengue antibody detection, we have evaluated dengue antigen from four individual sources. Those were NMRI, WRAIR, Rutgers University, and a Hawaiian biotechnology company (Hawaii Biotechnology Group, Inc.).

**(Figure 12)**

#### **Dengue Virus test**



The antigen sources varied in the following manner:

1. The Navy antigen (NMRI) was a native preparation partially purified by centrifugation. Most antigens processed in this manner do not work well in our type of format unless the typical antibody titers are very high. This preparation was not purified and contained media and animal sera straight from production. In this type of membrane format we normally require very pure specific proteins since we are observing direct binding and cannot use enzyme or other systems where the signal can be amplified. We saw very low reactivity with this preparation.
2. The Army (WRAIR) antigens are recombinant antigens for dengue 1 and 2. We do not know what part of the target has been produced (although we assume it is the entire envelope protein) but we do not see any reactivity with the clinical samples made available to us by WRAIR.
3. The antigens obtained from Rutgers University were small synthetic peptides of approximately 8 amino acid units long conjugated to BSA. I believe that



these were too short to be practical especially since you normally need about 15 amino acid residues to determine an epitope site and the chances of finding each antigenic site which patients may have titer to are slim. We have not seen any reactivity with these preparations in trials in our lab.

4. The final reagents are from a vaccine company in Aiea, Hawaii, Hawaii Biotechnology Groups, Inc., and are again recombinant viral antigens prepared in baculovirus. We did our initial testing with antigen from dengue-2, although we now have antigens for dengue 1, 3, & 4 which we have not as yet formatted on our assays. Dengue 1 and 2 are similar in nature so I am not surprised that we have already seen reactions on positive dengue 1 & 2 sera with the dengue 2 antigen. This antigen did produce good signals with all positive samples available to us from both WRAIR and the CDC. **(Information on those sera samples is provided as Attachment D).** Using a cruder preparation of the Hawaiian antigen did not work as well for sensitivity and we lost some of the positive samples. I believe we will need to stick with the purity level they initially supplied us with. Currently this assay picks up both IgG and IgM samples and we have shown 100% sensitivity and 100% specificity with 20 positive samples from the CDC and 2 serum pools. We have tested the assay using our panel of ten problematic negative samples and ten normal negative samples and have not observed any false positive results. We have sent samples to the CDC for further evaluation.

**Specific methodology for the dengue assay** - Using the antigen obtained from the Hawaii company, we applied at 1 mg/ml at 1 microliter per cm. The cards were blocked as the hanta assay. The same assay buffer was applied to the sample pads along with the same conjugate and amount of conjugate. We were not required to add any E coli extract or any other type of material such as detergents to the serum sample. However, in order to increase the assay time to five minutes using viscous samples, we found that it was necessary to dilute the sample 50:50 with passive gold diluent. If time for running and completing the test were not an important parameter, then neat serum was found to work just as well for specificity and sensitivity. The assay took up to 10 to 15 minutes for migration to occur with some samples when buffer was not added.

We may also, once the work is complete on all four dengue serotypes, format an assay with four separate bands representing those serotypes. Further development would aim to specifically identify separately the IgM and IgG reactions.

## **CONCLUSION**

We have shown that we can achieve equivalent visual detection using fluorescent markers, although we had originally thought that we could exceed the visual detection already standard in our colloidal gold format. We have concluded that the problem here lies within the limits of human visual acuity and that to detect the greater levels of

sensitivity that we still believe can be, or have been, achieved, that a specific UV reader is necessary. This would be fully investigated under a Phase II program. We have little doubt that we can eventually find the proposed 100-1000 fold increase in detection capability, bringing these simple assays to greater detection levels than current confirmatory testing capability.

Our reader system has met our expectations in this first go around. We have examined different readers and now feel confident that we can use this knowledge to modify existing, or engineer new, equipment that will result in a hand-held reader that can enhance interpretation of detection and provide quantification of results. In terms of ease of use, forward deployment, point of care, level of necessary expertise, this potential advancement in assay technology represents a clear advantage over the current body of knowledge and available equipment.

Temporarily discarding our efforts to approach a radial design for the multi-analyte format has enabled us to adjust our strip assay to the multi function design. We have not, as yet, worked enough with this assay to make any conclusive statements, but our previous experience with all of the component reagents bodes well for the eventual outcome.

Finally, the advancement of the hanta and dengue antibody detection assays has been a welcome product of our efforts. The early trials of both of these tests present us with a convincing potential for commercial marketability as well as DOD application.

## 3rd standard curve




Chart1

std100a




Chart1

hcg50a





Chart1

hcg0a

Actual scan  
of strip



— Series1

# **Arista Biologicals, Inc.**

115 Research Drive Bethlehem, PA 18015 phone(610)861-6901 fax(610)861-8247

ATTACHMENT B

March 17, 1997

Cheryl Lowery  
Contract Specialist  
R&D Contracts Branch B (CAL/dlh)  
USAMRAA  
Fort Detrick, MD 21702

Subject: SBIR Contract No. DAMD17-97-C-7027

Dear Ms Lowery:

We are interested in amending the above referenced SBIR contract to allow our company to receive and work with some additional unique military-relvant agents and reagents. WE can accomplish this revision to the original contract provisions within the current negotiated amount of \$98,150.00 and the Phase I completion date of May 24, 1997. The requested revisions are as follows:

- a. Delete the effort for enterotoxigenic *Escherichia coli* detection in stool specimens because the required reagents are not currently available.
- b. Add VEE virus detection. IN addition to usefulness for multi-analyte detection, this would be a good model for fluorescent labeling of the detector antibody.
- c. Add VEE virus IgM tests as a candidate for multi-test for IgM antibodies to dengue viruses and hantaviruses.
- d. Add plague F1 antigen detection as a candidate for multi-test for BW agents in serum.
- e. Add anthrax-protective antigen (PA) detection as a candidate for multi-test for BW agents in serum.
- f. Add ricin detection as a candidate for multi-test for BW agents in serum.
- g. Add Staphylococcal enterotoxin B (SEB) detection as a candidate for multi-test for BW agents in serum.

Thank you for your consideration in this matter. If you need more information, please contact me at 610-861-6901 or contact Dr. Greenspan at 908-730-8717 or fax to him at 908-735-6842.

Respectfully,

Dr. John Colanduoni  
President  
Arista Biologicals, Inc.

cc: Lt. Col. Roger Parker, USAMRIID



# Majesco DENGUEVIRUS

## A One Step Denguevirus Antibody Detection Test

*A rapid, one step test of the qualitative detection of IgM and IgG antibodies to Denguevirus in serum or plasma.*

*For professional in vitro diagnostic use only.*

### INTENDED USE

The Majesco One Step Denguevirus Test is a rapid chromatographic immunoassay for the qualitative detection of antibodies to Denguevirus in serum or plasma.

### SUMMARY

Denguevirus is a member of the flavivirus family, endemic in tropical and subtropical areas, and a major cause of morbidity and mortality in those regions. There are 4 identified sero-types of the virus (D1, D2, D3, D4), although the vast majority of infections are caused by just 2 (D1, D2), with many individuals exhibiting primary infection by either of the 2 major types and secondary infections by the other. The clinical spectrum of the disease ranges from asymptomatic to severe hemorrhagic/shock syndrome. Key to lowering morbidity and mortality is aggressive treatment in the early stages, with diagnosis dependent upon early identification of antibodies to the Denguevirus.

The Majesco Denguevirus One Step Test is a rapid test utilizing proprietary recombinant antigens to qualitatively detect the presence of both IgM and IgG antibodies to the D1 and D2 sero-types of the Denguevirus.

### PRINCIPLE

The Majesco Denguevirus One Step Test is a qualitative, solid phase, immunoassay for the detection of antibodies specific to Denguevirus 1 & 2 in serum or plasma. The membrane is pre-coated with Denguevirus recombinant antigens on the test line region and anti-conjugate on the control line region. During the testing, the serum or plasma sample reacts with the dye conjugate (anti-human antibody-colloidal gold conjugate) which has been pre-coated within the test device. The mixture then migrates upward on the membrane chromatographically by capillary action to react with the recombinant antigen on the membrane and generate a red line. Presence of this red line indicates a positive result, while its absence indicates a negative result. Regardless of the presence of detectable antibodies, as the mixture continues to migrate across the membrane to the immobilized anti-conjugate area, a red line at the control line region will always appear. The presence of this red line serves as verification for sufficient sample volume and proper flow and as a control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature (20-30°C). The kit is stable within the expiration date printed on the label. DO NOT FREEZE or use beyond the expiration date.

### MATERIALS PROVIDED

- 40 Majesco Denguevirus one step test devices
- 40 Disposable sample droppers
- 1 set of instructions for use

### MATERIALS REQUIRED BUT NOT PROVIDED

- Vacutainer tubes for serum or plasma collection
- Centrifuge

### ATTACHMENT C

- For professional in vitro diagnostic use only. Do not use after the expiration date
- Do not eat, drink, smoke in the area where the specimens and kits are handled
- Handle all specimens as though they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when assaying samples.

### SPECIMEN COLLECTION

- The Majesco Denguevirus One Step Test can be performed on serum or plasma only
- Remove the serum or plasma from the clot of red cells as soon as possible to avoid hemolysis. Only clear, non-hemolyzed specimens should be used.
- Testing should be performed as soon as possible after sample collection. Do Not leave samples at room temperature for prolonged periods. Specimens can be refrigerated at 2-8°C for up to 3 days. Otherwise specimens should be stored at below -20°C.
- Bring specimens to room temperature prior to testing. The frozen specimens must be completely thawed prior to testing. Specimens should not be repeatedly frozen and thawed.
- If specimens are to be shipped, they should be packed in compliance with Federal regulations covering the transportation of biologic agents.

### TEST PROCEDURE

Read the entire procedure carefully prior to performing any tests. **Allow test devices and serum samples to equilibrate to room temperature (20-30°C) prior to testing.**

- 1) Remove the Majesco Denguevirus test device from foil pouch (bring the test to room temperature before opening the pouch). Use device as soon as possible but within 1 hour after removal from pouch especially if the room temperature is more than 30°C and in high humidity environment.

2) Place the test device on a clean and level surface. Holding the dropper vertically, dispense three full drops of serum (~150ul) without air bubbles into the sample well "S" of the test device.

3) Wait for red lines to appear. The test should be read in approximately 10 minutes. It is significant that the background is clear before reading the test, especially when samples have low antibody titers, and only a weak line appears in the test region (T). Do not interpret results after 30 minutes.

### QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered as internal positive procedural control. A clear background in the result window is considered as internal negative procedural control. It is recommended that appropriate positive and negative serum or plasma sample controls be included in each days testing to verify proper test performance.

### INTERPRETATION OF RESULTS

**POSITIVE:** Two distinct lines will appear, one in the test region (T) and another in the control region (C).

**NEGATIVE:** Only a single red line appears in the control region (C). No apparent red or pink line appears in the test region (T).

**INVALID:** Control line fails to appear which means improper testing procedure or deterioration of reagents probably has occurred. In any event, repeat the test. If the problem persists, discontinue use of lot. Immediately and contact your local distributor.

**NOTES:** The shade of red color in the test line region (T) will vary depending on the concentration of antibodies present. However, neither the quantitative value nor the rate of increase in antibodies can be determined by this qualitative test.

### LIMITATION OF THE PROCEDURE

Majesco Denguevirus test is for in vitro diagnostic use only. This test should be used for the detection of antibodies to Denguevirus sero-types 1&2 in serum or plasma samples.

The test will only indicate the presence of antibodies to the specified Denguevirus strains in the specimen and should not be used as the sole criteria for the diagnosis of denguevirus infection. As with all diagnostic tests, a result must be considered in conjunction with other clinical information available to the physician. The Majesco Denguevirus test cannot detect extremely low concentrations of antibodies in specimens. If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is required. A negative result at any time does not preclude the possibility of Denguevirus infection.

### EXPECTED VALUES

Majesco Denguevirus test was compared to leading commercial EIA tests for denguevirus and the correlation between the systems was 98%.

### PERFORMANCE CHARACTERISTICS

#### Sensitivity

Majesco Denguevirus test was tested with a sensitivity panel containing sera positive to the D1 and D2 sero-types of the virus with resultant appropriate sensitivity to both strains.

#### Specificity

Recombinant antigens used for the Majesco Denguevirus test were developed in baculovirus. The two primary sero-types (D1, D2) produced positive results on the Majesco test. Specificity was tested using multiple well characterized samples of D1 & D2 with no incidence of false positive results.

### Precision

#### Intra-Assay

Within run precision was determined by using 15 replicates of three specimens (2 containing antibodies to D1&D2 and 1 without antibodies to either). The positive and negative values were correctly identified 98% of the time.

#### Inter-Assay

Between run precision was determined by using the same three specimens in 15 independent assays and with three different lots of Majesco Denguevirus test reaction cards over a 3 month period. The negative and positive values were correct 98% of the time.

### Reference Denguevirus Method

	positive	negative
<b>Majesco</b>	positive	145
	negative	5
<b>Method</b>	positive	0
	negative	150

### BIBLIOGRAPHY

1. Baron, EJ, Chang, RS, et al. Medical Microbiology. J. Wiley. 1994.

Manufactured for:

**Majesco Biologicals, Inc.**

115 Research Drive

Bethlehem, PA 18015

Tel: (610)861-6901

Fax: (610)861-8247

## ATTACHMENT D

### Sera Samples used to evaluate Dengue Assay

**NB:** All samples reacted with the assays using the Hawaii Bio. antigen.

From the DOD (WRAIR):

<u>Sample</u>	<u>Titre IgM</u>	<u>Titre IgG</u>
1041	50	600
1046	12800	1600
1068	12800	6400
1021	50	6400
1136	12800	6400
1146	50	400
1080	50	12800
1029	1600	1600
1138	12800	1600
1137	12800	12800
1122	50	1600
1037	50	12800
1123	100	400
1047	12800	1600

From the CDC:

**NB:** The samples from the CDC were sent in 2 shipments, 2 months apart, and were run in separate in-lab trials with different batches of tests and different batches of antigen reagent. Although there are some repeats in sample ID numbers, this amounts to 2 separate studies and the results, as given above, show consistency over time, from batch to batch of tests, and from batch to batch of reagent.

The heading designations in the first sampling group are defined as:

DPO = days post onset of symptoms

O.D. = optical density in the standard MAC-ELISA

First grouping of samples:

<u>Lab #</u>	<u>DPO</u>	<u>O. D.</u>	<u>Type of Sample</u>
334123	14	1.19	Primary Dengue 1
333870	5	0.72	Secondary Dengue 1
335445	10	0.26	"
335533	7	0.94	"
335704	17	0.93	"
294987	21	0.43	Primary Dengue 2
332823	8	0.35	Secondary Dengue 2
332960	8	0.53	"
332980	12	0.71	"
333797	11	1.83	"
335581	13	1.36	"
335927	16	0.56	"
332977	6	0.30	Undefined sera, D-1,2,4
332978	6	0.61	"
332983	5	0.65	"
333157	5	0.40	"
333163	6	0.39	"
333166	6	0.45	"
333286	5	0.31	"
333320	8	0.93	"
333353	7	0.29	"
333354	4	0.70	"
333358	8	0.76	"
333359	4	2.28	"

Second grouping of samples (ID#s only provided by CDC):

primary dengue 1

275558  
326250  
308224  
327327  
327483

secondary dengue 1

335545  
335533  
335704

primary dengue 2

334926  
336725  
294989

secondary dengue 2

334022  
335581  
335927  
336102  
336230

## ATTACHMENT E

### PREPARATION OF IMMUNOCHROMATOGRAPHIC STRIPS

#### TABLE OF CONTENTS

- 1) LIST OF MATERIALS SUPPLIERS
- 2) SOP 04 001.....PREPARATION OF BLOCK BUFFER
- 3) SOP 04 002.....PREPARATION OF SOLUTION "C" FOR THE  
TREATMENT OF ABSORBANT PADS
- 4) SOP 04 003.....PREPARATION OF SOLUTION "B" AND  
TREATMENT OF POLYESTER MEMBRANE FOR  
SPRAYING WITH GOLD CONJUGATE
- 5) SOP 04 004.....PREPARATION OF ANTIBODY SOLUTION FOR  
SPRAYING
- 6) SOP 04 005..... LAMINATION OF GLUE AND MEMBRANE TO THE  
PLASTIC CARD
- 7) SOP 04 006..... APPLICATION OF ASSAY REAGENTS TO THE  
MEMBRANE
- 8) SOP 04 007..... SPRAYING GOLD CONJUGATE ONTO POLYESTER  
2033
- 9) SOP 04 008..... ASSEMBLY OF THE DIAGNOSTIC TEST CARD
- 10) EXAMPLES OF BATCH RECORDS FOR MANUFACTURING

## LIST OF SUPPLIERS

### ABSORBANT PAPER

Schleicher and Schuell  
10 Optical Ave.  
P.O. Box 2012  
Keene N.H. 03431-2062  
FAX 603-357-3627  
TELEPHONE 800-245-4024  
Contact ; Nancy Moran

### CARDS (PLASTIC)

Estok Plastics co.  
434 Whitehead Rd.  
Trenton, N.J. 08619  
FAX 609-584-1833  
PHONE 609-584-1623  
609-586-4336

G&L  
Precision Die Cutting  
1766 Junction Ave.  
San Jose, CA 95112  
FAX 408-451-1199  
PHONE 800-327-4553  
408-453-9400

### CHEMICALS

Fisher Scientific  
FAX 800-926-1166  
PHONE 800-766-7000

Sigma Chemical Co.  
St. Louis, MO. 63178-9916  
FAX 800-325-5052  
PHONE 800-325-3010

### DESSICANTS

Grace-Davison Chemical Co.  
Curtis Bay, M.D.  
PHONE 410-354-8954  
Contact; Amy Lacher

DESSICANTS (continued)

United Dessicants  
101 Christine Drive  
Belen, N.M. 87002  
FAX 505-864-9296  
PHONE 505-864-6691

( we use  $\frac{1}{2}$  gm. packets)

Drierite  
W.A. Hammond Drierite Co. Ltd.  
Xenia, Ohio 45385  
937-376-2927

(this is the granular type )

GLUE

Adhesives Research Inc.  
P.O. Box 100  
Glen Rock, PA 17327  
FAX 717-235-8320  
PHONE 717-235-7979  
Contact; Karen Olsen

LAMINATOR

Bio.Dot. Inc. (see sprayers for address)

MEMBRANE

Millipore  
80 Ashby Road  
Bedford, MA 01730  
FAX 617-275-5550  
PHONE 800-225-3384  
Contact; William Mc Kensie, Marketing Manager  
Sheila, ext. 8662, Technical Service

POLYESTER 2033

Ahlstrom Filtration  
Scientific Specialties Group  
122 W. Butler St.  
P.O. Box A  
Mt. Holly Springs, PA 17065  
FAX 717-486-4878  
PHONE 800-326-1888  
Contact; Jack Ramsey

POUCHES

Acon Biotech (Hang Zhou) Co., LTD.  
148 Longjing Road  
Hang Zhou, ZheJiang  
P.R. of China 310013  
contact; Feng Lin  
( suppliers of large foil pouches)

POUCHES (continued)

Process Medical Systems  
25 East Loop Rd.  
Stonybrook, N.Y. 11790  
FAX 516-444-8845  
PHONE 516-444-8840

POUCHER

Bio.Dot., Inc. (see sprayers for address)

SLITTER/CUTTER

Bio. Dot., Inc. (see sprayers for address)

Kinematic Automation  
Twain Harte, CA 95383  
FAX 209-532-0248  
PHONE 209-532-3200  
contact; Ted Meigs

SPRAYERS

Bio. Dot., Inc. ( Bio Jet sprayer)  
17781 Sky Park Circle  
Irvine, CA 92714  
FAX 714 440- 3694  
PHONE 714-440-3685  
contact; John Geneau

Imagene Technology (Isoflow sprayer)  
P.O. Box 667  
Hanover, N.H. 03755  
FAX 603-298-8513  
PHONE 603 298-5632  
contact; Tim Chow



## ARISTA BIOLOGICALS

### ADDITIONAL POUCH MANUFACTURERS

- |  |  |
|--|--|
| <p>1. <b>ARVEY CORPORATION</b><br/>3500 North Kimball Avenue<br/>Chicago, IL 60618<br/>Tel No: (312) 463- 1400</p> <p>2. <b>FASER INDUSTRIES</b><br/>Medical Packaging Group<br/>Building P-1<br/>Andrea Boulevard<br/>Saddle Brook, NJ 07622<br/>Tel No: (201) 797-2600</p> <p>3. <b>KAPAK CORPORATION</b><br/>5305 Parkdale Drive<br/>St. Louis Park, MN 55416<br/>Tel No: (612) 541-0730</p> <p>4. <b>KENPAK</b><br/>5416 Union Pacific Avenue<br/>City of Commerce, CA 90022<br/>Tel No: (312) 463- 1400</p> <p>5. <b>LAMINATED FILM &amp;<br/>PACKAGING</b><br/>5070 Industry Drive<br/>Melbourne, FL 32940<br/>Tel No: (305) 259-2131</p> <p>6. <b>LPS INDUSTRIES</b><br/>113 North 13th Street<br/>Newark, NJ 07107<br/>Tel No: (201) 485- 4400</p> <p>7. <b>MODERN PACKAGING</b><br/>Shreve Street<br/>Mount Holly, NJ 08060<br/>Tel No: (609) 267-5900</p> <p>8. <b>RICHMOND<br/>TECHNOLOGY, INC</b><br/>1897 Colton Avenue<br/>Rodlands, CA 92374<br/>Tel No: (714) 794-2111</p> | <p>9. <b>RILEY &amp; GEEHR,<br/>INC.</b><br/>2205 Lee Street<br/>Evanston, IL 60202<br/>Tel No: (708) 869-8100</p> <p>10. <b>RJH &amp; ASSOCIATES, INC.</b><br/>419 Eisenhower Lane South<br/>Lombard, IL 60148<br/>Contact: Rock A. Hauptman<br/>Tel No: (708) 629-4010</p> <p>11. <b>SPEC-FAB COMPANY</b><br/>1818 Rowland Street<br/>Riverton, NJ 08077<br/>Tel No: (609) 829-5100</p> <p>12. <b>T.W KUTTER COMPANY</b><br/>91 Wales Avenue<br/>Avon, MA 02322<br/>Tel No: (617) 288-2600</p> |
|--|--|

ARISTA BIOLOGICALS. INC.

STANDARD OPERATING PROCEDURE

PAGE 1 OF 2

SOP 04 001 PREPARATION OF BLOCK BUFFER

```
=====
VERSION NO: 1.0                      SUPERSEDES: NEW FORMAT
PREPARED BY: B.SZILAGYI              DATE: 02/23/93
APPROVED BY: J.COLANDUONI            DATE: 02/23/93
REVIEWED BY: B.SZILAGYI              DATE: 09/04/95
=====
```

1.0 PURPOSE

1.1 To describe the preparation of block solution.

2.0 SAFETY

2.1 Refer to the MSDS sheet appropriate to the material being handled. No safety precautions are required.

3.0 EQUIPMENT

- 3.1 Analytical balance.
- 3.2 Glass beaker of adequate volume to contain buffer.
- 3.3 Magnetic stir plate.
- 3.4 Magnetic stir bar.
- 3.5 Graduated cylinder.
- 3.6 pH meter.

4.0 MATERIALS:

- 4.1 Albumin, bovine, Sigma A-7906, or equivalent.
- 4.2 Sucrose, Sigma S-8501, or equiv.
- 4.3 Phosphate buffered saline (PBS)  
[.01 M Na<sub>2</sub>HPO<sub>4</sub> dibasic, .15 M NaCl, pH 7.4]

5.0 PROCEDURE

- 5.1 Determine the volume of block solution needed.  
Note: 1000 ml. minimum for up to 50 strips membrane.
- 5.2 Pour entire volume of PBS into a glass beaker. Volume is same as final volume.
- 5.3 Weigh amount of albumin (BSA) to make 2.0% of final volume. Sprinkle BSA evenly over PBS and allow to dissolve without stirring.
- 5.4 Weigh amount of sucrose to equal 20% of final volume and stir until dissolved.
- 5.5 Adjust pH to 7.4 as necessary.
- 5.6 Store at 4-8 C refrigeration until use.

6.0 CLEAN UP

- 6.1 Wipe balance with damp cloth as necessary.
- 6.2 Clean used glassware by standard procedures.

ARISTA BIOLOGICALS, INC.

STANDARD OPERATING PROCEDURE

PAGE 1 OF 2

SOP 04 002      PREPARATION OF SOLUTION "C" FOR TREATMENT  
OF ABSORBANT PADS

```
=====
VERSION NO:  2.0                      SUPERSEDES: 6/24/93
PREPARED BY: B.SZILAGYI                DATE:  8/19/93
APPROVED BY: J.COLANDUONI              DATE:  8/19/93
REVIEWED BY: B.SZILAGYI                DATE:  9/07/95
=====
```

#### 1.0 PURPOSE

1.1 This procedure provides guidelines for the preparation of SOLUTION C and the treatment of absorbant paper.

#### 2.0 SAFETY

2.1 Protective eyewear and clothing must be worn as indicated by the MSDS of the chemicals used.

#### 3.0 EQUIPMENT

- 3.1 Analytical balance
- 3.2 pH meter
- 3.3 Magnetic stir plate and stir bar
- 3.4 Beaker
- 3.5 Tray

#### 4.0 MATERIALS

- 4.1 Triton X-100, Sigma X-100 or equivalent
- 4.2 Sodium Tetraborate, Fisher S249-500, or equiv.
- 4.3 HCl 4N or NaOH to adjust pH
- 4.4 DI water
- 4.5 Absorbant paper 8"x10", Grade 903, Schleicher & Schuell
- 4.6 Gloves

#### 5.0 PROCEDURE

- 5.1 Weigh the desired amount of 0.1 M SODIUM BORATE:  
0.1 x 381.37 (formula wt) x \_\_\_\_\_liters final  
vol.=  
\_\_\_\_\_grams sodium borate. Lot  
#\_\_\_\_\_.
- 5.2 Dissolve the above in final volume less 100 mls DI water.
- 5.3 Add the desired amount of TRITON X-100 10%: Final

volume in mls x 0.10 = vol. of Triton X-100 to be added.

\_\_\_\_\_ mls. x 0.10 = \_\_\_\_\_ mls. Triton X-100.

Lot # \_\_\_\_\_.

- 5.4 Dissolve thoroughly with mixing.
- 5.5 Adjust the pH to 10.0.
- 5.6 Pour solution into a tray and briefly but thoroughly soak the absorbant paper sheet one at a time in the solution.  
Wear gloves when handling the sheets.
- 5.7 Remove the sheet from the solution and blot the sheet between two dry absorbant sheets.
- 5.8 Airdry the sheets completely on racks, taking care not to deform the sheets.
- 5.9 Cut the sheets into strips 3/4" x 180 mm.

#### 6.0 CLEAN UP

- 6.1 Store unused SOLUTION C in capped container at room temperature.
- 6.2 Rinse tray with DI water.
- 6.3 Store cut strips in a clean sealed container.

## STANDARD OPERATING PROCEDURE

SOP 04 003

```
=====
VERSION NO: 2.0                                SUPERSEDES: 06/24/93
PREPARED BY: B.SZILAGYI                        DATE: 08/19/93
APPROVED BY: J.COLANDUONI                      DATE: 08/19/93
REVIEWED BY: B.SZILAGYI                       DATE: 09/07/95
=====
```

1.1 The purpose of this procedure is to provide guidelines for the preparation of SOLUTION B and the treatment of polyester membrane grade 2033 (Ahlstrom).

2.1 Protective eyewear and clothing is to be worn as suggested by the MSDS of the chemical indicated.

- 3.1 Analytical balance
- 3.2 pH meter
- 3.3 Magnetic stir plate and stir bar
- 3.4 Beaker
- 3.5 Tray
- 3.6 Gloves

- 4.1 Polyvinyl Alcohol (PVA) 30,000 - 70,000 MW, Sigma P8136
- 4.2 Sodium Phosphate, dibasic, Sigma S-0876, or equiv.
- 4.3 Triton X-100, Sigma X-100
- 4.4 Bovine Serum Albumin (BSA), Sigma A-7906, or equiv.
- 4.5 DI water
- 4.6 Polyester membrane, grade 2033, Ahlstrom Filtration

5.1 Weigh the desired amount of POLYVINYL ALCOHOL  
0.5%.  
Dissolve thoroughly with the application of

moderate heat with stirring. Lot # \_\_\_\_\_.

Final volume in mls.  $\times 0.005 =$  \_\_\_\_\_ grams PVA.

5.2 Weigh the desired amount of SODIUM PHOSPHATE 50 mM.

$142 \times .001 \times 50 \times \text{final vol. in liters} = \text{gm Na}_2\text{HPO}_4$

7.1  $\times$  \_\_\_\_\_ L = \_\_\_\_\_ gm  $\text{Na}_2\text{HPO}_4$ . Lot # \_\_\_\_\_

Remove from heat and dissolve with mixing.

5.3 Add the desired amount of Triton X-100 1.0%.

Final volume in mls  $\times 0.01 =$  mls Triton X-100.

\_\_\_\_\_ mls  $\times 0.01 =$  \_\_\_\_\_ mls Triton. Lot # \_\_\_\_\_

Mix thoroughly.

5.4 Cool solution to less than 40C.

5.5 Weigh the desired amount of BSA 5.0%.

Final volume in mls.  $\times 0.05 =$  gms BSA.

\_\_\_\_\_ mls  $\times 0.05 =$  \_\_\_\_\_ gms BSA. Lot # \_\_\_\_\_

5.6 Add BSA to cooled solution and mix gently until dissolved.

5.7 Adjust pH to 7.4.

5.8 Pour Solution B into a tray and, wearing gloves, dip the polyester sheets into the solution to wet completely. Soak 1 hour at room temperature.

5.9 Shake off excess solution and blot the sheets between paper towels. Dry on racks.

5.10 Store treated sheets in a clean, covered container until ready to spray with the gold conjugate.

#### 6.0 CLEAN UP

6.1 Unused Solution B may be stored under refrigeration at 4-8 C for up to three days.

SOP 04 004 PREPARATION OF ANTIBODY SOLUTION FOR SPRAYING

=====

VERSION NO: 1.0

SUPERSEDES: NEW FORMAT

PREPARED BY: B.SZILAGYI

DATE: 02/23/93

APPROVED BY: J.COLANDUONI

DATE: 02/23/93

REVIEWED BY: B.SZILAGYI

DATE: 09/09/95

=====

1.0 PURPOSE

- 1.1 To prepare the antibody solution for spraying onto membrane.

2.0 SAFETY

- 2.1 Wear protective clothing, gloves, and goggles as appropriate to the nature of the antibody.

3.0 EQUIPMENT

- 3.1 Glass pipets and bulb  
3.2 Graduated cylinder

4.0 MATERIALS

- 4.1 Phosphate Buffered Saline (PBS)  
(0.01 M Na<sub>2</sub>HPO<sub>4</sub>, dibasic, + 0.15 M NaCl, pH 7.4)  
4.2 Stock antibody

5.0 PROCEDURE

- 5.1 Calculate the final volume of antibody solution needed:

5.1.1 
$$\text{Volume of antibody solution in mls.} = \frac{\text{number of cards to be sprayed} \times \text{microliters of antibody per card}}{\text{divided by 1000.}}$$

- 5.2 Dilute the stock antibody to desired concentration (AC) in PBS:

5.2.1 
$$\text{AC divided by the concentration of stock in mg/ml} \times \text{final volume of antibody solution in mls.} = \text{vol. of stock antibody in mls. needed.}$$

5.2.2 
$$\text{Volume of PBS in mls.} = \text{final volume in mls.} - \text{vol. of stock antibody in mls.}$$

6.0 CLEAN UP



- 6.1 Wash cylinder in nonfilming laboratory soap, deproteinize between different antibody types with 0.1 NaOH,, and rinse thoroughly in DI water.
- 6.2 Place stock antibody, and the diluted working concentration of antibody solution under refrigeration at 4 - 8 C.

SOP 04 005 LAMINATION OF GLUE AND MEMBRANE TO PLASTIC CARD

=====

VERSION NO: 2.0

SUPERSEDES: 2/23/93

PREPARED BY: B.SZILAGYI

DATE: 09/09/95

APPROVED BY: J.COLANDUONI

DATE: 09/09/95

REVIEWED BY:

DATE:

=====

1.0 PURPOSE

- 1.1 To describe the lamination of the glue and membrane components to the plastic card.

2.0 SAFETY

- 2.1 Observe proper safety precautions while using the cutter.

3.0 EQUIPMENT

- 3.1 Paper cutter or equivalent  
3.2 Template

4.0 MATERIALS

- 4.1 White polystyrene cards, 2 1/4" x 180 mm, Estok Plastics.  
4.2 Arcare 7841 1" double-coated 1 mil polyester film, liner both sides, Adhesives Research.  
4.3 Arcare 7841 1/2" double-coated 1 mil polyester film, liner both sides, Adhesives Research  
4.4 Arcare 7841 3/4" double-coated 1 mil polyester film, liner both sides, Adhesives Research  
4.5 Membrane, High Flow Nitrocellulose, 1" x 180mm, unspecified pore size, unbacked, Millipore

5.0 PROCEDURE

- 5.1 Cut the Arcare rolls into strips 180 mm long.  
5.2 Remove the protective liner from one side of the 3/4" strip.  
5.3 Using the template, place the tape, glue side down, on the plastic card along the bottom edge of the card.  
5.4 Remove the liner from one side of the 1" tape and place it on the card 3/4" from the bottom edge, taking care not to overlap with the 3/4" tape.

- 5.5 Remove the protective liner from one side of the 1/2" tape and place it glue side down along the top edge of the card just above the 1" tape.
- 5.6 Cut strips of membrane 1" x 180 mm long.
- 5.7 Wearing gloves to protect the membrane from oils, remove the top liner from the 1" tape and position the membrane onto the glue. Take care in positioning the membrane because once in place, it cannot be removed without damaging the membrane. Make sure there are no bubbles in the membrane and avoid dust particles.

NOTE: Do not rub the surface of the membrane during positioning; instead cover the positioned membrane with a second plastic card and apply even pressure.

#### 6.0 CLEAN UP

- 6.1 Place all unused tape into plastic bags for storage.
- 6.2 Store unused membrane in a protective container.
- 6.3 Store finished cards in a dust-free container until use.

SOP 04 006 APPLICATION OF ASSAY REAGENTS TO MEMBRANE

```
=====
VERSION NO: 1.0                      SUPERSEDES: NEW FORMAT
PREPARED BY: B.SZILAGYI              DATE: 09/09/95
APPROVED BY: J.COLANDUONI            DATE: 09/11/95
REVIEWED BY:                        DATE:
=====
```

1.0 PURPOSE

- 1.1 To describe the application of antibody/antigen reagents to membrane, blocking and drying procedures.

2.0 SAFETY

- 2.1 Gloves should be worn at all times to protect both the technician and the membrane.
- 2.2 Protective clothing should be worn.

3.0 EQUIPMENT

- 3.1 Sprayer, Biodot, or equivalent
- 3.2 Incubator/dessicator 37C
- 3.3 Tray for blocking

4.0 MATERIALS

- 4.1 Prepared working concentration of antibody/antigen solution.
- 4.2 Cards, 2 1/4" x 180 mm, with laminated membrane
- 4.3 Paper towels
- 4.4 Blocking solution
- 4.5 Gloves
- 4.6 Compressed gas supply

5.0 PROCEDURE

- 5.1 Program the appropriate setting onto the sprayer.
- 5.2 Place the card on the platform and align the spray head over the desired line of spray. Secure the card with magnets.
- 5.3 Fill the resevoir with the antibody/antigen solution to be sprayed.
- 5.4 Turn on the gas supply and begin the spray, watching the spray pattern closely for consistency and accuracy.
- 5.5 Dry the sprayed cards on a rack in the incubator for 30 minutes at 37C.

- 5.6 Soak the cards face down in block solution for 30 minutes at 37C. Wear gloves when handling cards in block solution.
- 5.7 Remove the cards from the block solution and blot between paper towels membrane side down, taking care not to rub the membrane. Insufficient blotting may affect function of test strips.
- 5.8 Dry the blocked cards overnight at 37C then store in desiccated containers at room temperature. Optimal relative humidity in storage container should be 20% or less.

#### 6.0 CLEAN UP

- 6.1 Thoroughly clean the spraying syringes and jets as recommended by the manufacturer.
- 6.2 Place all unused reagents into refrigeration at 4 - 8C.
- 6.3 Daily discard all unused block solution.

ARISTA BIOLOGICALS, INC.

STANDARD OPERATING PROCEDURE

PAGE 1 OF 2

SOP 04 007      SPRAYING GOLD CONJUGATE ONTO POLYESTER 2033

```
=====
VERSION NO: 1.0                      SUPERSEDES: NEW FORMAT
PREPARED BY: B.SZILAGYI              DATE: 08/19/93
APPROVED BY: J.COLANDUONI            DATE: 08/19/93
REVIEWED BY: B.SZILAGYI              DATE: 09/09/95
=====
```

1.0 PURPOSE

- 1.1 To describe the process for spraying gold conjugate.

2.0 SAFETY

- 2.1 Gloves and protective clothing should be worn.

3.0 EQUIPMENT

- 3.1 Sprayer, Biodot or equivalent
- 3.2 Incubator/dessicator set at 37C

4.0 MATERIALS

- 4.1 Prepared gold conjugate suspension
- 4.2 Polyester sheets, 2033 grade, treated with Solution C.
- 4.3 Gloves
- 4.4 Compressed gas supply

5.0 PROCEDURE

- 5.1 Program the appropriate settings on the sprayer.
- 5.2 Check alignment of spray head with desired line of spray.
- 5.3 Secure polyester sheet on the platform with magnets.
- 5.4 Turn on the gas supply and begin spray, watching spray pattern closely for consistency and accuracy.
- 5.5 Sprayed lines of gold conjugate may be repeated at 1/4" intervals until the sheet is full.
- 5.6 Dry thoroughly at 37C no longer than overnight.
- 5.7 Cut into 1/4" strips such that the line of gold conjugate run along the top of each strip.

## 6.0 CLEAN UP

- 6.1 Thoroughly rinse the syringe with water, then clean with 1% Triton X-100 (Camag) or 0.02% Bioterge (BioDot), then rinse completely with DI water. Use approximately 10 volumes each of the detergent and the final water rinse.
- 6.2 Wipe sprayer platform with a damp cloth to remove conjugate residues.

ARISTA BIOLOGICALS, INC.

STANDARD OPERATING PROCEDURE

PAGE 1 OF 2

SOP 04 008

ASSEMBLY OF DIAGNOSTIC TEST CARDS

```
=====
VERSION NO:  1.0                      SUPERSEDES: NEW FORMAT
PREPARED BY:  B.SZILAGYI              DATE: 02/23/93
APPROVED BY:  J.COLANDUONI            DATE: 02/23/93
REVIEWED BY:  B.SZILAGYI              DATE: 09/09/95
=====
```

1.0 PURPOSE

- 1.1 To provide guidelines for the assembly of diagnostic test cards.

2.0 SAFETY

- 2.1 Not applicable.

3.0 EQUIPMENT

- 3.1 Template

4.0 MATERIALS

- 4.1 Cards with laminated membrane, sprayed and blocked  
4.2 5/8" x 180 mm strips of treated absorbant paper  
4.3 3/4" x 180 mm strips of treated absorbant paper  
4.4 1/4" strips of sprayed gold conjugate  
4.5 gloves

5.0 PROCEDURE

NOTE: The following assembly procedures must be performed under conditions of low humidity, ideally less than 40%. Desiccate top and bottom pads for two hours at 37 C before assembly.

- 5.1 Remove the liner from the 3/4" glue strip along the bottom of the card. Wear gloves to apply all components.
- 5.2 Position the gold strip with the line of gold conjugate along the top such that there is approximately a 1 - 1.5 mm overlap onto the membrane.
- 5.3 Place the 3/4" absorbant pad along the bottom edge of the card, taking care to leave the gold strip exposed.
- 5.4 Remove the liner from the 1/2" glue strip along the top of the card. Place the 5/8" absorbant strip along the top of the card overlapping onto the membrane by 1 to 1.5 mm.



- 5.5 Store the finished cards dessicated at room temperature until ready for cutting into test strips. Relative humidity of storage container must be 20% or less.
- 5.6 If the strips are to be placed into cassettes and pouched, remove only the number of strips from the container that can be placed in fifteen minutes.

CLEAN UP

- 6.1 Place all unused materials in their original containers.

Name of Assay: \_\_\_\_\_ Lot #: \_\_\_\_\_ Number of Cards: \_\_\_\_\_

<b>CONTROL LINE:</b> Reagent: _____ Lot #: _____ Stock Concentration: _____ Working Concentration: _____ Dilution made by: _____ Date _____ Stock Solution: _____ Equilibration: _____ PBS: _____ Batch Volume: _____ Sprayed by: _____ Date _____ Dispenser: _____ Dispense Volume: _____	<b>TEST LINE:</b> Reagent: _____ Lot #: _____ Stock Concentration: _____ Working Concentration: _____ Dilution made by: _____ Date _____ Stock Solution: _____ Equilibration: _____ PBS: _____ Batch Volume: _____ Sprayed by: _____ Date _____ Dispenser: _____ Dispense Volume: _____	<b>GOLD LINE:</b> Reagent: _____ Lot #: _____ Stock OD <sub>540</sub> : _____ Working OD <sub>540</sub> : _____ Stabilized By: _____ Date _____ Stock Solution: _____ Equilibration: _____ PGB: _____ Sprayed by: _____ Date _____ Dispenser: _____ Dispense Volume: _____
--	---	---

<b>TOP PAD:</b> Lot #: _____ Cut by: _____ Size: _____	<b>MEMBRANE:</b> Lot #: _____ Laminated by: _____ Block Buffer made by: _____ Treated by: _____ Date _____	<b>POLYESTER SHEET:</b> Lot #: _____ Solution B made by: _____ Treated by: _____ Date _____ Cut by: _____	<b>BOTTOM PAD:</b> Lot #: _____ Solution C made by: _____ PH = _____ Treated by: _____ Date _____ Cut by: _____ Size: _____
---	--	---	---


CARD ASSEMBLY by: \_\_\_\_\_ Date \_\_\_\_\_  
STRIP CUTTING by: \_\_\_\_\_ Date \_\_\_\_\_  
DEVICE ASSEMBLY by: \_\_\_\_\_ Date \_\_\_\_\_  
DESICCANT Lot #: \_\_\_\_\_  
DROPPER LOT #: \_\_\_\_\_  
POUCHING by: \_\_\_\_\_ Date \_\_\_\_\_

TOTAL NUMBERS OF TESTS: \_\_\_\_\_  
QUALITY CONTROL by: \_\_\_\_\_ DATE: \_\_\_\_\_

# Membrane Placement

Place Membrane Here

# Chromatographic Assay Card Template



— Antibody Capture Line (Printed on Nitrocellulose Membrane)

# Membrane Placement


	<b>Place Membrane Here</b>	

# Top Absorbant Pad Cutting Template

## Bottom Absorbant Pad Cutting Template

# Absorbant Pad Placement

Align Top Pad To Line



Align Bottom Pad To Line



5.1 Location \_\_\_\_\_  
5.2 Room temperature \_\_\_\_\_

ARISTA BIOLOGICALS, INC.  
STANDARD OPERATING PROCEDURE AND BATCH RECORD

SOP 06 007      PREPARATION OF NORMAL PHOSPHATE SALINE  
FROM 10 X PBS

=====

VERSION 1.0	SUPERSEDES: NEW FORMAT
PREPARED BY: B.SZILAGYI	DATE: 10/16/95
APPROVED BY:	DATE:

=====

1.0 PURPOSE

1.1 To provide guidelines and documentation for the preparation of normal phosphate buffered saline from 10X PBS.

2.0 SAFETY

2.1 Not applicable

3.0 BATCH VOLUME \_\_\_\_\_ PREPDATE \_\_\_\_\_ TECH \_\_\_\_\_

MATERIALS	MANUFACTURER & LOT#	STD AMT (L)	ACT AMT
10X PBS		100 mls./ L	
DI water		900 mls.	
NaOH 4N		as necessary	

4.0 PROCEDURE

- 4.1 Based on batch size calculate the amount of 10X PBS needed. Measure the calculated amount in a graduated cylinder and pour into an adequately sized flask.
- 4.2 Calculate the volume of DI water needed and pour that measured amount into the flask.
- 4.3 Stir contents until thoroughly mixed, approximately twenty minutes.

5.0 STORAGE LOCATION AND CONDITIONS:

- 5.1 Location \_\_\_\_\_
- 5.2 Store covered at room temperature for up to two weeks. Discard if there is evidence of microbial contamination.



DEPARTMENT OF THE ARMY  
US ARMY MEDICAL RESEARCH AND MATERIEL COMMAND  
504 SCOTT STREET  
FORT DETRICK, MARYLAND 21702-5012

REPLY TO  
ATTENTION OF:

MCMR-RMI-S (70-1y)

4 Dec 02

MEMORANDUM FOR Administrator, Defense Technical Information  
Center (DTIC-OCA), 8725 John J. Kingman Road, Fort Belvoir,  
VA 22060-6218


SUBJECT: Request Change in Distribution Statement

1. The U.S. Army Medical Research and Materiel Command has reexamined the need for the limitation assigned to technical reports written for this Command. Request the limited distribution statement for the enclosed accession numbers be changed to "Approved for public release; distribution unlimited." These reports should be released to the National Technical Information Service.

2. Point of contact for this request is Ms. Kristin Morrow at DSN 343-7327 or by e-mail at Kristin.Morrow@det.amedd.army.mil.

FOR THE COMMANDER:

Encl

  
PHYLLIS M. RINEHART  
Deputy Chief of Staff for  
Information Management

ADB218773	ADB229914
ADB223531	ADB229497
ADB230017	ADB230947
ADB223528	ADB282209
ADB231930	ADB270846
ADB226038	ADB282266
ADB224296	ADB262442
ADB228898	ADB256670
ADB216077	
ADB218568	
ADB216713	
ADB216627	
ADB215717	
ADB218709	
ADB216942	
ADB216071	
ADB215736	
ADB216715	
ADB215485	
ADB215487	
ADB220304	
ADB215719	
ADB216072	
ADB222892	
ADB215914	
ADB222994	
ADB216066	
ADB217309	
ADB216726	
ADB216947	
ADB227451	
ADB229334	
ADB228982	
ADB227216	
ADB224877	
ADB224876	
ADB227768	
ADB228161	
ADB229442	
ADB230946	
ADB230047	
ADB225895	
ADB229467	
ADB224342	
ADB230950	
ADB227185	
ADB231856	